

Product Datasheet

CD59 Antibody (aa29-102) (orb1539403)

Description	CD59 Antibody (aa29-102)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	A recombinant human CD59 protein 29-102 aa (BC001506).
Target	CD59
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide. PBS, pH 7.4, 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA (1:20000 - 1:80000), IHC, IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.