

Product Datasheet

PRICKLE2 Antibody (orb1539411)

Description

PRICKLE2 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P, WB
Immunogen	A 15 AA synthetic peptide near the center of human PRICKLE2. The immunogen is located within aa 420-470 of PRICKLE2.
Target	PRICKLE2
Preservatives	PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C for up 3 months. Store at -20°C for up to 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: PRICKLE2 antibody can be used for detection of PRICKLE2 by Western blot at 1 ug/mL. Antibody can also be used for immunofluorescence starting at 20 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in mouse samples and Immunofluorescence in rat samples. All other applications and species not yet tested. Western Blot:
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IF (20 µg/ml), IHC, IHC-P (5 µg/ml), WB (1 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558