

Product Datasheet

**MCUR1 / CCDC90A Antibody (aa252-301)
(orb1539433)**

Description

MCUR1 / CCDC90A Antibody (aa252-

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide located between aa252-301 of human MCUR1 (Q96AQ8, NP_001026883). Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Marmoset (100%); Galago, Mouse, Rat, Guinea pig, Platypus (92%); Elephant, Dog, Horse, Turkey, Chicken (85%).
Target	MCUR1 / CCDC90A
Preservatives	PBS, 0.09% sodium azide, 2% sucrose. PBS, 0.09% sodium azide, 2% sucrose
Concentration	0.5 mg/ml
Storage	Short term: Store at 2-8°C. Long term: Aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Western Blot: Suggested dilution at 1µg/ml in 5% skim milk / PBS buffer, and HRP conjugated anti-Rabbit IgG should be diluted in 1: 50,000 - 100,000 as second antibody. ELISA titer in peptide based assay: 1:1562500.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), WB (0.2 - 1 µg/ml)
Expiration Date	12 months from date of receipt.