

Product Datasheet

ITGA2B / CD41 Antibody (aa122-171) (orb1539437)

Description	ITGA2B / CD41 Antibody (aa122-171)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide located between aa122-171 of human ITGA2B (P08514, NP_000410). Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Monkey, Elephant, Dog (100%); Galago, Marmoset, Horse, Pig (92%); Guinea pig (85%); Mouse, Rat, Rabbit (84%).
Target	ITGA2B / CD41
Preservatives	PBS, 0.09% sodium azide, 2% sucrose. PBS, 0.09% sodium azide, 2% sucrose
Concentration	0.5 mg/ml
Storage	Short term: Store at 2-8°C. Long term: Aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), WB (1 µg/ml)
Expiration Date	12 months from date of receipt.