

## Product Datasheet

**KRT18 / CK18 / Cytokeratin 18 Antibody (clone C-04, Biotin) (orb1539479)**

**Description**

KRT18 / CK18 / Cytokeratin 18

|                            |                                                                                                                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                |
| <b>Reactivity</b>          | Mammal                                                                                                                                                                                               |
| <b>Conjugation</b>         | Biotin                                                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA, FC, ICC, IHC, IHC-P, IP, WB                                                                                                                                                                   |
| <b>Immunogen</b>           | Cytoskeleton preparation of epidermal carcinoma cell line A431                                                                                                                                       |
| <b>Target</b>              | KRT18 / CK18 / Cytokeratin 18                                                                                                                                                                        |
| <b>Preservatives</b>       | PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide                                                                                                                                     |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                                                                              |
| <b>Storage</b>             | Store at 2°C to 8°C. Do not freeze.                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                                |
| <b>Application notes</b>   | Further information: Suggested working dilution is 1:500. Indicated dilution is recommended starting point for use of this product. Working concentrations should be determined by the investigator. |
| <b>Isotype</b>             | IgG1                                                                                                                                                                                                 |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                           |
| <b>Dilution Range</b>      | ELISA, Flo (1:500), ICC, IHC, IHC-P (10 µg/ml), IP, WB                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                      |