

Product Datasheet

CD68 Antibody (clone KP1) (orb1539490)

Description	CD68 Antibody (clone KP1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IHC-P, IP, WB
Immunogen	Subcellular fraction of human alveolar macrophages.
Target	CD68
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C or at -20°C. Store undiluted. Avoid freeze-thaw cycles. Microcentrifugation recommended if solution contains precipitate.
Note	For research use only
Application notes	Further information: Flow Cytometry: Use 10ul of the suggested working dilution to label 106 cells in 100ul.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo, IF, IHC, IHC-P (5 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.