

Product Datasheet

GSK-3 Beta mouse Monoclonal Antibody HRP Conjugated (orb1539567)

Description	GSK-3 Beta mouse Monoclonal Antibody HRP Conjugated
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ICC, IHC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	P49841
Dilution Range	WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200
Expiration Date	12 months from date of receipt.