

Product Datasheet

Bax mouse Monoclonal Antibody PE-Cy7 Conjugated (orb1539570)

Description	Bax mouse Monoclonal Antibody PE-Cy7 Conjugated
Species/Host	Mouse
Conjugation	PE/Cy7
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	Q07812
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.