

Product Datasheet

Myc tag mouse Monoclonal Antibody AP Conjugated (orb1539608)

Description	Myc tag mouse Monoclonal Antibody AP Conjugated
Species/Host	Mouse
Conjugation	AP
Tested Applications	ICC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	P01106
Dilution Range	WB=1:10000-20000 ICC=1:50-200
Expiration Date	12 months from date of receipt.