

Product Datasheet

Myc tag mouse Monoclonal Antibody Bio Conjugated (orb1539609)

Description	Myc tag mouse Monoclonal Antibody Bio Conjugated
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ICC, IF, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	P01106
Dilution Range	WB=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.