

Product Datasheet

phospho-Smad2/Smad3 (Thr8) Polyclonal Antibody PE Conjugated (orb1539683)

Description	phospho-Smad2/Smad3 (Thr8) Polyclonal Antibody PE Conjugated
Species/Host	Rabbit
Conjugation	PE
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P84022
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50- 200 IF=1:50-200
Expiration Date	12 months from date of receipt.