


Product Datasheet

ITGA6 antibody (orb153986)

Description	Rabbit polyclonal antibody to ITGA6. Integrins are alpha / beta heterodimeric cell surface receptors...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Integrin alpha 6 (ITGA6/CD49f)(NP_001073286.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125KDa
Uniprot ID	P23229
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...