

## Product Datasheet

# HSPG2/Heparan Sulfate Proteoglycan 2 Polyclonal Antibody FITC Conjugated (orb1539937)

|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | HSPG2/Heparan Sulfate<br>Proteoglycan 2 Polyclonal<br>Antibody FITC Conjugated  |
| <b>Species/Host</b>        | Rabbit                                                                          |
| <b>Conjugation</b>         | FITC                                                                            |
| <b>Tested Applications</b> | IF                                                                              |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                           |
| <b>Clonality</b>           | Recombinant                                                                     |
| <b>Source</b>              | Rabbit                                                                          |
| <b>Uniprot ID</b>          | <b>P98160</b>                                                                   |
| <b>Dilution Range</b>      | IF=1:50-200                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |