

Product Datasheet

CD133 Polyclonal Antibody FITC Conjugated (orb1539990)

Description	CD133 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	O43490
Dilution Range	Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.