

Product Datasheet

Cytokeratin 19 Polyclonal Antibody PE-Cy7 Conjugated (orb1540424)

Description	Cytokeratin 19 Polyclonal Antibody PE-Cy7 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy7
Tested Applications	FC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P08727
Dilution Range	Flow-Cyt=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.