

Product Datasheet

Placental alkaline phosphatase (PLAP) Polyclonal Antibody FITC Conjugated (orb1540570)

Description	Placental alkaline phosphatase (PLAP) Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P05187
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.