

Product Datasheet

SAA1/Serum Amyloid A Polyclonal Antibody FITC Conjugated [Out of stock] (orb1540598)

Description	SAA1/Serum Amyloid A Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P0DJ18
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.