

Product Datasheet

CK14/17/42/10 Polyclonal Antibody FITC Conjugated (orb1540675)

Description	CK14/17/42/10 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P02533
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.