

Product Datasheet

Phospho-Estrogen Receptor alpha (Tyr537) Polyclonal Antibody Bio Conjugated (orb1540744)

Description	Phospho-Estrogen Receptor alpha (Tyr537) Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	IF, IHC-Fr, IHC-P, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P03372
Dilution Range	B=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.