

Product Datasheet

GAD65 / GAD67 Polyclonal Antibody Bio Conjugated (orb1540967)

Description	GAD65 / GAD67 Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ICC, IF, IHC
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q05329
Dilution Range	IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.