


Product Datasheet

Occludin antibody (orb154101)

Description	Rabbit polyclonal antibody to Occludin. This gene encodes an integral membrane protein which is loca...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Occludin(NP_002529.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	Q16625
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...