

Product Datasheet

HLA-DR/HLA DRB1 Polyclonal Antibody PE-Cy5.5 Conjugated (orb1541071)

Description	HLA-DR/HLA DRB1 Polyclonal Antibody PE-Cy5.5 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P01903
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.