

Product Datasheet

Nuclear Pore Complex Proteins/Nup107 Polyclonal Antibody Bio Conjugated (orb1541100)

Description	Nuclear Pore Complex Proteins/Nup107 Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ICC, IF, IHC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P57740
Dilution Range	WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.