

Product Datasheet

Aquaporin 5 Polyclonal Antibody FITC Conjugated (orb1541170)

Description	Aquaporin 5 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P55064
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.