

Product Datasheet

Collagen I antibody (FITC) (orb15412)

Description

Collagen I antibody (FITC)

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Human, Mouse, Rat, Sheep

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

Synthetic peptide (please contact us for details)

Target

Collagen I

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

161 kDa

Uniprot ID
P02458, P11087, P02454
NCBI
NM_001844.4, NP_001835.3
Entrez
1280
Dilution Range

IF/ICC:1:50-200 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.