

Product Datasheet

HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated (orb1541211)

Description	HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy5
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P30443
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.