


Product Datasheet

PLAUR antibody (orb154129)

Description	Rabbit polyclonal antibody to PLAUR. Urokinase plasminogen activator surface receptor (uPAR, also kn...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PLAUR(NP_002650.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36-72kDa
Uniprot ID	Q03405
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...