

Product Datasheet

FoxP3 Polyclonal Antibody APC Conjugated (orb1541301)

Description	FoxP3 Polyclonal Antibody APC Conjugated
Species/Host	Rabbit
Conjugation	APC
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9BZS1
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.