

Product Datasheet

alpha smooth muscle Actin Polyclonal Antibody PE Conjugated (orb1541307)

Description	alpha smooth muscle Actin Polyclonal Antibody PE Conjugated
Species/Host	Rabbit
Conjugation	PE
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P62736
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50- 200 IF=1:50-200
Expiration Date	12 months from date of receipt.