

Product Datasheet

Myc tag Polyclonal Antibody HRP Conjugated (orb1541396)

Description	Myc tag Polyclonal Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-Fr, IHC-P, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P01106
Dilution Range	WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200
Expiration Date	12 months from date of receipt.