

Product Datasheet

Eukaryotic aspartyl protease Polyclonal Antibody PE Conjugated (orb1541484)

Description	Eukaryotic aspartyl protease Polyclonal Antibody PE Conjugated
Species/Host	Rabbit
Conjugation	PE
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	G0YYY2
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.