

Product Datasheet

CD133 Polyclonal Antibody FITC Conjugated (orb1541591)

Description	CD133 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC
Immunogen	751-848 AA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Purity	Affinity purified by Protein A
Source	Rabbit
Uniprot ID	O43490
Dilution Range	Flow-Cyt=1:50-200
Expiration Date	12 months from date of receipt.