

Product Datasheet

GFAP Polyclonal Antibody Bio Conjugated (orb1541635)

Description	GFAP Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ICC, IF, IHC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P14136
Dilution Range	WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.