

Product Datasheet

AIF mouse Monoclonal Antibody FITC Conjugated (orb1541699)

Description	AIF mouse Monoclonal Antibody FITC Conjugated
Species/Host	Mouse
Conjugation	FITC
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	095831
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.