

Product Datasheet

Nestin Polyclonal Antibody RBITC Conjugated (orb1541711)

Description	Nestin Polyclonal Antibody RBITC Conjugated
Species/Host	Rabbit
Conjugation	RBITC
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P21263
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.