

Product Datasheet

SIRT6(Phospho-Ser338) Antibody HRP Conjugated (orb1542111)

Description	SIRT6(Phospho-Ser338) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q8N6T7
Dilution Range	Western blotting: 1:100-1000Immunohistochemistry1:100-500
Expiration Date	12 months from date of receipt.