

Product Datasheet

TrkA (Phospho-Tyr496) + TrkB (Phospho-Tyr516) + TrKC (Phospho-Tyr516) Antibody FITC Conjugated (orb1542115)

Description	TrkA (Phospho-Tyr496) + TrkB (Phospho-Tyr516) + TrKC (Phospho-Tyr516) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q16620
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.