


Product Datasheet

smo antibody (orb154214)

Description	Rabbit polyclonal antibody to smo. The protein encoded by this gene is a G proteincoupled receptor t...
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SMO(NP_005622.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86KDa
Uniprot ID	Q99835
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brai...