


Product Datasheet

SOS1 antibody (orb154220)

Description	Rabbit polyclonal antibody to SOS1. Son of sevenless (SOS) was first identified in Drosophila as a g...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SOS1(NP_005624.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	175kDa/200kDa
Uniprot ID	Q07889
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...