

Product Datasheet

PLC gamma 2 (Phospho-Tyr759) Antibody FITC Conjugated (orb1542403)

Description	PLC gamma 2 (Phospho-Tyr759) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P16885
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.