

Product Datasheet

PERK(Phospho-Thr980) Antibody HRP Conjugated (orb1542410)

Description	PERK(Phospho-Thr980) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9Z2B5
Dilution Range	Western blotting: 1:100-1000Immunohistochemistry1:100-500
Expiration Date	12 months from date of receipt.