

Product Datasheet

5-Lipoxygenase(Phospho-Ser663) Antibody FITC Conjugated (orb1542417)

Description	5-Lipoxygenase(Phospho-Ser663) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P09917
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.