

Product Datasheet

GSK3 Alpha Beta(Phospho-Tyr279+Tyr216) Antibody FITC Conjugated (orb1542726)

Description	GSK3 Alpha Beta(Phospho-Tyr279+Tyr216) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P49840
Dilution Range	FCM: 1:20- 100Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.