

Product Datasheet

PD-1 CD279 Antibody FITC Conjugated (orb1542791)

Description	PD-1 CD279 Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q15116
Dilution Range	Immunofluorescence: 1:50-200 FC=1:50-200
Expiration Date	12 months from date of receipt.