

Product Datasheet

Integrin Alpha V + Beta 3 Antibody Cy5 Conjugated (orb1543035)

Description	Integrin Alpha V + Beta 3 Antibody Cy5 Conjugated
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P05106
Dilution Range	FCM:1:20-100IF:1:50-200
Expiration Date	12 months from date of receipt.