

## Product Datasheet

### TAZ Antibody HRP Conjugated (orb1543094)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | TAZ Antibody HRP Conjugated                               |
| <b>Species/Host</b>        | Rabbit                                                    |
| <b>Conjugation</b>         | HRP                                                       |
| <b>Tested Applications</b> | IHC-P, WB                                                 |
| <b>Storage</b>             | Store at 4°C for 12 months.                               |
| <b>Note</b>                | For research use only                                     |
| <b>Clonality</b>           | Recombinant                                               |
| <b>Source</b>              | Rabbit                                                    |
| <b>Uniprot ID</b>          | <b>Q9GZV5</b>                                             |
| <b>Dilution Range</b>      | Western blotting: 1:100-1000Immunohistochemistry1:100-500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |