



Product Datasheet

YAP antibody (orb154318)

Description	Rabbit polyclonal antibody to YAP. The transcriptional coactivator Yes-associated protein (YAP) has ...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human YAP1(NP_001123617.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72kDa
Uniprot ID	P46937
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...



Immunohistochemical staining of human ga...