

## Product Datasheet

### ROR1 Antibody FITC Conjugated (orb1543431)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | ROR1 Antibody FITC Conjugated   |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Conjugation</b>         | FITC                            |
| <b>Tested Applications</b> | IF                              |
| <b>Storage</b>             | Store at 4°C for 12 months.     |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Recombinant                     |
| <b>Source</b>              | Rabbit                          |
| <b>Uniprot ID</b>          | <b>Q01973</b>                   |
| <b>Dilution Range</b>      | Immunofluorescence:<br>1:50-200 |
| <b>Expiration Date</b>     | 12 months from date of receipt. |