

Product Datasheet

AKT1 2 3 (Phospho-Ser472 Ser473 Ser474) Antibody HRP Conjugated (orb1543489)

Description	AKT1 2 3 (Phospho-Ser472 Ser473 Ser474) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P47196
Dilution Range	Western blotting: 1:100-1000Immunohistochemistry1:100-500
Expiration Date	12 months from date of receipt.