


Product Datasheet

TNF alpha antibody (orb154376)

Description	Mouse monoclonal antibody to TNF alpha
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant human TNF-alpha
Target	TNF alpha
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Functional application: Neutralization. ELISA: It can be used as capture antibody in combination with biotinylated antibody MAb11.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P01375
Entrez	7124
Dilution Range	Functional application: Neutralization. ELISA: It can be used as capture antibody in combination with biotinylated antibody MAb11.
Expiration Date	12 months from date of receipt.



ELISA analysis of TNF alpha (0 - 4.1 ng/...